## Multi-variable indicators

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## Results Reference Guide

<p>Turns to  or lighter</p>	←		<b>Unprocessed</b>
<p>Turns to  or lighter</p>	←		<b>Failures</b>
<p>Turns to  or lighter</p>	←		<p>One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again. The main failure causes are:</p> <ul style="list-style-type: none"> <li>• Inadequate cycle duration</li> <li>• Inadequate cycle selection</li> <li>• Overloading / incorrect loading of the sterilizer</li> <li>• Improper sterilization agent concentration</li> <li>• Sterilizer malfunction</li> <li>• Blocking indicator with large metal mass</li> </ul>
<p>Turns to  or lighter</p>	←		<b>Processed and Satisfactory</b>
<p>Turns to  or lighter</p>	←		

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**ATTENTION:** If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through *Bionova*® Software. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).

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**NOTE:** Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

## Multi-variable Indicators

For Plasma or Vaporized Hydrogen Peroxide sterilization

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Type 4 Indicators according to ISO 11140-1:2014

### Product description

Chemdye® CD40 Multi-variable Indicators have been designed to react to Plasma or Vaporized Hydrogen Peroxide sterilization processes within loads, ensuring an adequate control of the effectiveness of sterilization processes. Purple ink was developed to change to green when the process reaches the stated values of the critical process variables. Multi-variable indicators allow to rapidly detect whether the critical parameters of the sterilization process for which they have been designed to respond have been reached.

### Characteristics

Chemical strips (105 x 18 mm) printed with indicator ink.  Metals free

### Indications for use

Use Chemdye® CD40 Multi-variable Indicators for monitoring Plasma or Vaporized Hydrogen Peroxide sterilization processes.

Conditions: 6 minutes, 50 °C, 2,3 mg/l H<sub>2</sub>O<sub>2</sub>.

### Precautions

Do not expose Chemdye® CD40 Multi-variable Indicators to Steam, Ethylene Oxide or any sterilization process other than Plasma or Vaporized Hydrogen Peroxide.

### Instructions for use

- Place a Chemdye® CD40 Multi-variable Indicator within each package, bag or tray that will be sterilized by Plasma or Vaporized Hydrogen Peroxide.
- In rigid containers place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- After completing the process, analyze results. The indicator must turn to the reference color for considering that the indicated conditions were met.

### Storage and shelf life

- Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.
- After use, the indicator will not change within 1 year when stored at previously indicated conditions.

### Disposal

Discard chemical indicators after use according to your country's healthcare and safety regulations.

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